


# SURGICAL DISPLAYS



# iM-MS270P

## 27"

FULL HD, AR ANTI-GLARE  
PROTECTIVE FRONT GLASS,  
EXQUISITE APPEARANCE



[www.imed-inc.com](http://www.imed-inc.com)



MADE IN CHINA



# iM-MS270P



## 27"

### FULL HD, AR ANTI-GLARE PROTECTIVE FRONT GLASS, EXQUISITE APPEARANCE

The iM-MS270P is a 27-inch surgical monitor that features a 1920x1080P IPS HD screen, high-transmittance AR anti-glare protective front glass, and a mirror appearance design on the outside.

This design effectively prevents dust and liquid pollution during surgeries, meeting high-level purification requirements of operating rooms.

The product is equipped with multiple signal interfaces, making it suitable for various image displays. Its standard installation structure is compatible with operating room equipment.





**True Color  
Reproduction**



**LED Backlight**



**Exquisite and  
Real Image**



**Full View  
Displaying Panel**



**Comply with International  
Medical Standard**

## IPS

Display Size	27"
Color	1.07B
Resolution	1920×1080
Display Ratio	16 : 9
Brightness	1000 cd/m <sup>2</sup>
Contrast Ratio	1000 : 1 (TYP)
Viewing Angle	178°/178°
Response Time	14ms (TYP)
Input Ports	VGA/ DVI/ S-VIDEO/ VIDEO/ BNC(SDI)
Output Ports	DVI/ BNC(SDI)
Power Supply	DC 24V
Power Consumption	≤80W max
MTBF	>30,000h
Structural Appearance	White Metal Shell With Protective Front Glass
Dimension	647*409*63 mm (W*H*D) (Without Stand)
Weight	8.5±0.5 kg (Without Stand)
Installation Standard	VESA 100x100 mm
Application	Digital Video-endoscope System / ENT
Certifications	CE



# CONTACT US



## **Medical Monitors Industry; Made in China**

ADDRESS OF LINES PRODUCTION:

Longhua District, Shenzhen, Shenzhen, Guangdong  
Province, China,



Tel/Fax: +971 43337118

Sales-Team: +971 58 582 4100

Technical Support-Team: +971 58 546 4100



[www.imed-inc.com](http://www.imed-inc.com)



[iMedRadiology](#)



[imed-inc-company](#)

---