

MOBILE X-RAY SYSTEMS



iMMD-110




www.imed-inc.com

 **MADE IN TURKEY**



iMed.

Inspire Your Mind



iMMD-110

**High Quality Image
With Simple
Operation**



**Focal Film Distance
Range: 30 - 120 (CM)
120° Rotation Tube**



**Wide Range
Monoblock
Rotation**




The Soft keypad, Digital display of KV and mAs value on the operating console enable exposure parameters to set easily and quickly.
The anatomical programs (APR) features with symbolic buttons, Contribute to intuitive operation with a minimum number of key stroke

TECHNICAL SPECIFICATIONS

Output Power:	5.6 KW
Max mA:	100 mA
MAx Kv:	110 Kv
Input Power:	220v
Time Range:	0.01 Sec to 6 Sec
mAs Range:	0.4 mAs to 200 mAs
Control Method:	HF-Inverter Method

X-Ray Tube

Max. Anode Heat Capacity:	42.8 kHU(30kj)
X-Ray Tube Facal Spot Size:	1.8 mm
Target Angle:	15





CONTACT US



H. Rifat Pasa Yüzer Havuz Sk. Perpa Ticaret Merkez B Blok Kat:13
No:2552 34384 SiSli /İstanbul Türkiye



Tel/Fax: +971 43337118
Sales-Team: +971 58 582 4100
Technical Support-Team: +971 58 546 4100



www.imed-inc.com



[iMedRadiology](#)



[imed-inc-company](#)
